




Product - Information



Kaps Balancing System

Ideal Balance for your Microscope

The Kaps Balancing System allows a better and easier movement of the microscope.

Die 45° inclinable coupling is mounted between the microscope carrier and the swivel arm.

Through additional accessories as secondary observation or video camera arise asymmetrical load allocation. With the Kaps balancing system this is no longer a problem. With the help of the coupling you can optimally balance out the microscope head, so that it doesn't come to undesirable turning movement.

With a swing screw you can adjust the motion of the turning movement.

- Order no. 112-420-000 Balancing unit (for 101-130-000 / 101-132-000)
- Order no. 112-450-000 45° inclinable balancing coupling