



Product - Information



Kaps - photo tube

for digital (DSLR) cameras (Canon or Nikon)

The photo tube digital allows various DSLR cameras to be connected to Kaps SOM[®]-series, slit lamps, colposcopes as well as iris diagnostic microscopes

- Multicoated antireflection high end optics for high transmission ≥ 90
 - Full-size, distortion free digital camera image without vignetting
 - Easy and fast connection to the beam splitter „plug and play“ with height and side correct image orientation and individually upright adjustment, which allows easy and fast changing of the camera position by turning
 - Direct connection possibility to digital cameras via T2-adaptation without adjustment, thus nearly all DSLR-s are applicable
 - Only 320 grams, thus there is no impairment of the microscope handling
 - Through the compact design it does not produce undesirable turning movements
 - Attachable camera types: Canon and Nikon (Identify type of camera before ordering)
 - Adaptation possible to Zeiss (other brands on request)
 - Additional accessories on request
- ➡ Order no 101-320-100

