



Product - Information



Kaps - Free Motion magnetic brake system

The free motion magnetic brake provides easy, quick and accurate positioning of the microscope head.

The free motion magnetic brake in combination with the balancing equipment allows the microscope head to be positioned in an easy, effortless way. The procedure is simple. The two magnetic brakes are activated and released by simply pressing a button. The microscope head can now be moved effortlessly, precisely and safely into any working position using only one hand. Once the buttons are released, the position is fixed.

This also works when using a maximum of accessories, such as an adaptation of a co-observation device or a digital camera.

The Kaps SOM® series maintains the perfect balance through high mechanical quality, easy positioning and a unique price-performance ratio.

- ➔ order no 112.870-000 in conjunction with microscopes with a motor function
- ➔ order no 112.890-000 in conjunction with microscopes with manual functions