

Kaps HD Photo Adapter for SONY NEX-5/-6/-7 and Alpha 5000/6000

Discover the HD documentary with its unique variety of details

- Apochromatic optics latest generation, for brilliant and true color images and videos
- Highest resolution for current and future high-definition camera systems with APS-C-Format
- Ideally matched to camera and microscope-Full format without any annoying vignetting shading
- Compact overall height and minimum weight, no adverse effect on the word with the microscope
- Live full HD-video streaming in 1920 x 1080 p via HDMI
- Smooth connection of the camera to a PC-Documentation of the highest quality
- Adapter for connection to Zeiss microscopes is available (Order no 101-328-000)



