



Product - Information


Kaps - Media Work Station BE

System for the collecting and archiving of video findings in the medical field

Based on state-of-the-art computer technology
Compatible with Windows 7 / VISTA / XP / DV and USB digital Video

Multilingual - English, French, German, Italian, Russian, Spanish, Turkish, Chinese PR

Intuitive User interface with static layout, self-explanatory work flow
Mobile & Touch



MANAGEMENT

- ➔ Patient data, images and comments are stored in structured database
- ➔ Automatic assignment of patient ident numbers
- ➔ Search and sort functions for patient data
- ➔ Full text search function Automatic backup with scheduled interval

FINDINGS

- ➔ Multiple-split display during recording
- ➔ Recording controllable by foot or camera
- ➔ Video playback with a slide bar for individual selection of time positions
- ➔ Input of diagnoses comments

DOCUMENTATION

- ➔ Graphics marks for emphasizing of findings
- ➔ Automatic creation of customizable HTML reports with images
- ➔ Automatic cropping of black frames in printouts
- ➔ Direct print of images with embedded patient data and other additional information

HOTLINE SUPPORT

With Maintenance Agreement, reaction within 1 working day



Product - Information

Kaps - Media Work Station BE

Basic Package

ITEM	DESCRIPTION	ARTICLE #
Software	Software for HD/SD Video Documentation, storing, viewing and analysis of patient video findings, session, finding management	
Kaps-Media-Workstation BE Version 2.x	Full HD video, PAL/NTSC video Finding Archiving, photo editor, graphical labeling. HTML report with patient data and images, direct printing of patient data. User interface „touch optimized“ Compatible with Windows 8/7 and Vista Multilingual (DE/EN/FR/IT/ES/RU/TR) Medical device class 1	
Maintenance Agreement	One-Year Maintenance Agreement Telephone- and Email-support, remote maintenance Guaranteed response within 1 business day	
Installation Medium	Kaps - Media Work Station BE Pack Bag Installation data carrier in a self-sealing paper	
User Manual	Printed Version for the actual Software Version	
Package	Basic Package Kaps - Media Work Station BE Contains all the items above For initial order, a complete package must be ordered	800-150-000

Options

ITEM	DESCRIPTION	ARTICLE #
Maintenance Agreement (sold directly from Software Partner)	One-Year Maintenance Agreement Telephone- and Email-support, remote maintenance Guaranteed response with 1 bussiness day	800-150-001
Network License (only for networking several Kaps - Media Work Stations among them self, otherwise not required)	Network Software Option Central storage for up to 5 jobs, simultaneous video recording and play- back on multiple workstations Central user administration, supporting key License is required for each workstation in the network	800-150-002
Dicom - Push	Dicom - Push option Dicom export of medical reports and pictures via Dicom - Push	800-150-003
Foot - Switch	Foot - Switch Option Compatible to Windows 8/7/VISTA/SP without driver installation Controls image AND video recording of the Media Work Station software simultaneously	800-150-004

OUR FOCUS IS ON YOUR VISION