



Product - Information



Kaps - beam splitter with endoscope adapter

for retrofitting any diagnosis-microscopes with a endoscope camera

Kaps beam splitter with endoscope adapter is compact and efficient.

- very quick and easy to install - apply between enlargement changer and binocular tube - without adjustment - ready
 - only 34,3 mm height, therefore only a slight increase in length of microscope body
 - only 300 grams in weight, therefore it does not affect movement of the microscope body
 - high quality optics ensures an excellent depth of focus and contrast as an measurement
 - fits to all endoscope connectors
 - can be adapted for use on competition microscopes (Zeiss / Leica)
- ➔ Order no. 101-890-000