



## Product - Information



## Kaps - Vario Objective 200 - 400 mm

### Experience new Ergonomics at Workplace - with Kaps Vario Objective

The Kaps Vario Objective stands for ergonomics at workplace. The objective is simple to use and very comfortable. The optics allows you to move up and down the focal plane from 200 - 400 mm without changing lenses only by turning a knob.

The Vario optics design delivers extraordinary image quality, improved resolution and distortion free image with a high depth of focus. Precision „Made in Germany“ - this was the reason why Leica Camera manufactured the complete optics and coated it with one of the highest multi-layer antireflection coatings.

For all motorized operation microscope the vario objective is also available as a motorized version.

- Best ergonomics with all your treatments
- Premium apochromatic optics with multi coating
- Multi-layer antireflection coatings
- Stay focused quickly and easily, with a turn of the knob
- Focusing range  $f = 200-400$  mm

➡ Order no. 162-700-000 (manual version)

➡ Order no. 162-720-000 (motorized version)





## Product - Information

### Technical Data Vario Objective 200 - 400 mm

Binokular tube  $\phi$  -210° f=182 mm

Ocular WF 12,5 x

Zoom: 1:6

#### Zoom factors

Distance	0.3	0.63	1	1.6	2
----------	-----	------	---	-----	---

#### Magnification field of view mm

200	2.0	4.3	6.8	10.9	13.6
	71	37	22	15	12
250	1.8	3.7	5.9	9.5	11.8
	80	40	25	16	14
300	1.6	3.3	5.2	8.4	10.5
	92	46	28	18	16
350	1.4	3.0	4.7	7.5	9.4
	103	51	32	21	18
400	1.3	2.7	4.3	6.8	8.5
	125	58	36	23	20