

Kaps - 45° Ergonomic Wedge with Beam Splitter and Binocular Turntable

An unbeatable combination for a high level of convenience at the workplace.

Dentists, in particular, often complain about back or neck problems. Thanks to the Kaps 45° ergonomic wedge, including binocular turntable, these are now a thing of the past.

The user always sits comfortably in an upright, ergonomically correct position and can work in a much more relaxed way.

Even with varying working positions, or if the microscope is pivoted away to the side, the user always keeps the patient in focus. A simple hand movement and the preferred working position can be assumed again. The neck and back are noticeably relieved.



45° Ergonomics Wedge with Beam Splitter

The 45° ergonomics wedge with beam splitter is simply inserted between the magnifying changer/zoom and binocular turntable.



45° ergonomics wedge unilateral 50:50

Order no 101-990-100

45° ergonomics wedge unilateral 80:20

Order no 101-991-100



45° ergonomics wedge bilateral 50:50 / 80:20

Order no 101-993-100

Binocular Turntable

- The low height of just 20 mm ensures a minimal extension of the microscope head
- It weighs only 240 g. This means that no harmful interference occurs during handling of the microscope
- The rotation angle is $\pm 24^{\circ}$ and is rasterized in $4 \times 6^{\circ}$ steps
- With the clamping screw, locking of the position is possible in any position
- It is always possible to adapt all Kaps microscopes at a later point
- An adaptation of competitor models is also possible
- Order no. 101.989-000

